

QUANTUM

3D MICROSCOPE™

Powered by IBIS® Technology

Take firearm and tool mark identification to the next level

➤ Bring objectivity to common source conclusions



ULTRA | Forensic
Technology

QUANTUM

3D MICROSCOPE™

For Firearm and Tool Mark Identification



The Quantum 3D Microscope (Q3M) has capabilities beyond those of conventional comparison microscopes. Technological advancements in 3D topography measurements are unlocking much-needed innovations.

Q3M provides firearm and tool mark examiners with the best 3D visual and quantitative tools for common source determination.



Bullets and Other Small Objects

The Q3M S1 model captures the rifling marks on fired bullets, and tool marks on other small objects, such as press marks on pills and tablets, as well as chamber, extractor, and ejection port marks on cartridge cases.



Support Expert Conclusions Using Objective Methods

Firearm and tool mark identification is evolving and now requires 3D measurements to support expert conclusions using objective methods that provide confidence levels and error rates.

Powered by IBIS® Technology

Q3M capabilities are derived from 3D technology innovations in IBIS Search Networks that were developed over more than a decade.


While IBIS focuses on finding the needle-in-the-haystack (previously unrelated firearm crimes), Q3M focuses on specific casework and studies, evaluating and quantifying the strength of agreement in common source determinations.



Crime evidence



Suspect firearm recovered



Test-fired bullets




3D Acquisition



Q3M can save hours or days compared to doing the same casework on a comparison microscope.

Reach more conclusions with high-confidence in a fraction of the time.


3D Comparison Viewing



3D Shape View

3D Surface View

Quantitative Analysis




3D Comparison Viewing

A whole new level of viewing is made possible with 3D microscopy. More details can be seen with less effort.




Features

- Surface and shape views, always in focus
- Up to 6 objects compared simultaneously side-by-side
- Hairline moveable across overlapped surfaces
- Horizontal and vertical side-by-side views
- Overlap, move, and rotate in locked or free modes
- Sectional line profiles show precise topography
- Visual enhancements emphasize the marks
- Bullet best match assistance aligns the similarities
- Simulated lighting and surface reflectivity options
- Graphically annotate comparisons



Bullet comparisons



Ejection port mark comparison

Quantitative Analysis


Topography measurements feed algorithms that offer objective results in support of scientific studies and expert conclusions with confidence levels and error rates.

Compare bullets from unknown and known sources, and graphically represent the strength of common source similarity using Q3M's innovative RBL Graph¹.



RBL Graph

- The RBL method visually demonstrates quantifiable differences between matching and non-matching conditions.
- Correlation algorithms provide line counting and pattern matching scores (LCS and PMS).
- Best Phase mode is well-suited for intact bullets. Best LEA mode is well-suited for partial or damaged bullets.
- False Match Rate (FMR) provides a reliable error rate to support expert conclusions. Currently for conventional rifling, with plans to expand.
- Subclass Warning Coefficient indicates the possible influence of subclass characteristics.
- Clustering automatically groups same source items based on FMR as well as phase consistency (patent pending).




The False Match Rate (FMR) for a given similarity score represents the probability that two bullets that were not fired from the same firearm would generate a greater score.

¹ Roberge, D., Beauchamp, A., & Lévesque, S. (2019). Objective Identification of Bullets Based on 3D Pattern Matching and Line Counting Scores. *International Journal of Pattern Recognition and Artificial Intelligence*, 33(11)

3D Microscope


Technical Characteristics

- Non-Linear Photometric Stereo 3D sensor capable of measuring purely specular surfaces and diffusive surfaces (patented)
- Pilot cameras with smart positioning function
- Intelligent surface tracking for full object wraparound or across fragments
- Surface width: 2.8 mm
- Lateral resolution: 2.98 µm/pixel
Depth resolution: Less than 0.3 µm
- Calibration to traceable measurement standards
- Object dimensions: up to 50 mm long and 28 mm diameter, including all bullet calibers (0.17 to 0.50)




Customized Mounting Stubs

A set of 10 customized stubs is provided for the easy mounting of objects of various shapes and sizes. All calibers of bullets, in various shapes, can be mounted, including test fires, and damaged and fragmented evidence.



Stub is magnetically secured to shaft and automatically retracts

Set of 10 customized stubs




+ Add-On to IBIS BULLETTRAX™

The Q3M Acquisition Unit can be shared as part of an IBIS BULLETTRAX Acquisition Station for bullet entries onto an IBIS Search Network.

This can be advantageous if the use of the equipment can accommodate the volume of IBIS bullet entries and 3D microscope casework.

Flexible Data Management



- Organize captured data and results in workspaces
- Easily document the work product in notes
- Output images, screenshots, notes, and results for expert reports and external information management systems
- Import and export 3D acquisition data in native and X3P/OpenFMC formats



Standalone Workstation and Remote Extensions

Although the standard Q3M is a standalone workstation, acquisitions can be provided to remote workstations via workspace files.

For example, a technician can acquire the physical evidence from multiple cases so that examiners can work in parallel on comparison viewing and quantitative analysis. Examiners can share workspaces for peer reviews.



Acquisition Unit
Dimensions (W x H x D):
61.7 cm x 37.5 cm x 57.9 cm
(24.3" x 14.8" x 22.8")
Acquisition Unit Weight:
33 kg (72 lb.)

Training, Proficiency Testing, and Research Studies

- Acquire objects once, and freely view and share
- Build training sets and reference sets
- 3D topography data provides reliable measurements
- Perform studies to advance the firearm and tool mark identification discipline
- Collaborate in community research projects based on shared 3D data

QUANTUM

3D MICROSCOPE™

Powered by IBIS® Technology

Innovating
today
for a safer
tomorrow

ULTRA. | Forensic
Technology

Follow us



www.ultra-forensictechnology.com